


IPRO-Software

## winIPRO Interface with Zeiss i.Profiler®




Zeiss i.Profiler (EyeExam) 23.03.2009, 09:35

	Sph	Cyl	A°	R	L
R	+0,48	-0,25	158	8,05 8,49 e0,59	8,06 8,29 8,96
L	+0,63	-0,24	15	176° 8,18 7,60	e0,78 8,06 8,34

Beckham, C


**NEW:** i.Profiler data can now also be called up in the winIPRO Eye Examination Report. Here, the radii are shown additional to the basic refraction data.



### Simply a new tab in winIPRO

A simple click on the winIPRO tab i.Profiler leads directly to a new dimension. With this new interface, you can import the data measured with the i.Profiler easily into winIPRO.


The ordering of the i.Scription lenses at Zeiss is available at the touch of a button, and finally, you can manage the delivery values in winIPRO. As usual, a lot of know-how from the IPRO developers has been put into this, and of course IPRO has been involved in the developing process at Zeiss from the beginning, to find the best solution for interfacing the i.Profiler with other IT. Please order the new interface from your IPRO-consultant at the same time as you order the device at Zeiss, to be able to process lens ordering at Zeiss with the new measurement data via winIPRO.



ZEISS device i.Profiler

i.Profiler uses the so-called wavefront analysis to create an individual vision profile of the patient's eyes. During the test, real-life demands for the eye are considered and calculated, because life is not just black-and-white and fixed – it is eventful, with change of light, reflexions, colours, change of contrast etc.





#### Job procedure with the Zeiss i.Profiler

- ▶ Patient is entered or opened in winIPRO.
- ▶ Data from i.Profiler are retrieved Input or automatic transfer of refraction values into the spectacles order.
- ▶ Data transfer to other consultation or measurement programs, such as Winfrol, RVT, Video Infral etc, and storing of the results.
- ▶ If desired, i.Scription optimisation can be selected during lens selection.
- ▶ Lenses are ordered as usual with winIPRO.
- ▶ Enter the optimum values in delivery values tab

## IPRO is on the cutting edge of interface technology

If our partners from the industry develop new devices and software, we are present and develop the fitting interfaces in close co-operation. Sometimes this is new ground for us, too, but very often our developers' ideas make it into the design of these innovations. We are proud of these research investments, because our customers can rely on winIPRO to be on the cutting edge of technology. To find out more about IPRO-interfaces, please refer to the "Factsheet Interfaces" at [ipro.eu/modules/factsheets](http://ipro.eu/modules/factsheets)

#### winIPRO i.Profiler Interface

##### The Features:

Data retrieval from the i.Profiler®, storing in the spectacles order and transfer during lens ordering.  
Storing of the delivered values.

*winIPRO Ordering System and winIPRO Eye Examination Management are required*



For more Info about winIPRO Interface with Zeiss i.Profiler®, please contact your IPRO consultant.

Service and Support have a face

## This is ours - winIPRO

**Accept our invitation and contact us.  
All our team look forward to it.**

For more info contact:

Jay Greenslade  
j.greenslade@iprouk.com  
Phone 07595964146

German HQ:

info@ipro.de  
Phone: +49 (0) 7152 9333 0  
Fax: +49 (0) 7152 9333 30  
www.ipro.de

IPRO GmbH  
Steinbeisstr. 6  
71229 LEONBERG  
GERMANY

IPRO UK c/o Haag Streit UK Ltd  
Edinburgh Way, Harlow, Essex, CM20 2TT

[www.iprouk.com](http://www.iprouk.com)